





Product Specifications RW 24 V2 – 07/12/2016

PRONATIV® 95 INSTANT

Native Whey Protein Isolate

Description

PRONATIV® 95 Instant is a native whey protein isolate obtained by membrane filtration. Thanks to its instant quality PRONATIV® 95 is easy to rehydrate.

COMPOSITION (/ 100	g powder)	
Moisture (%)	max	6.0
Protein, as is (N*6.38) (%)	min	89.0
Protein, dry basis (N*6.38) (%)	min	94.0
Fat (%)	max	1.9
of which Lecithin added (%)	max	1.5
Lactose (%)	max	3.0
Minerals (%)	max	3.0

Applications

Dairy ingredient particularly adapted to sports nutrition, dietetics and clinical nutrition thanks to its amino acid composition (leucine)

Packaging

25 kg bags, multiwall kraft paper bag with inner polyethylene bag liner

Shelf Life and Storage

18 months in non-opened bag in a cool and dry place away from all sources of odours and tastes

Suggested Labelling

Milk proteins, emulsifier: soy lecithin

SENSORIAL
Colour: cream
Taste and Odour: neutral

CHEMICAL AND PHYSICAL			
Nitrates* (ppm)	max 50		
Phosphatase*	Negative		
Scorched particles ADPI (/25 g)	Disc A		
pH*	6.7 ± 0.3		
Bulk density*	0.35 ± 0.05 g/ml		
Wettability	max 20 s at 25°C		
* Indicative values			

BACTERIOLOGICAL (CFU*)			
Total Plate Count	/1g	< 10 000	
Yeasts & Moulds	/1g	< 50	
Coliforms	/1g	< 10	
Salmonella	/ 375 g	None	
Staphylococcus coag. + (including S. aureus)	/ 1 g	None	
Listeria	/ 25 g	None	
Contact us for other analysis Analysis according to standard methods (IDF or ISO) * Colony-forming unit			

NUTRITIONAL VALUES (calculated)			
Energy Value (/ 100 g) : 1600 kJ / 380 kcal			
Nutritional Information *			
Sugar (Lactose)	3.0 g		
Saturated Fatty Acids 0.4 g			
Sodium 100 mg			
*Theoretical values/100 g			

MINERALS (mg/100g powder)*		
Calcium	300	
Magnesium	53	
Phosphorus	200	
Potassium	400	
Iron	0.6	
Chlorides	180	
*Indicative values		

All rights to the information contained herewith belong to Lactalis Ingredients. The information is confidential and may not be disclosed to third parties or exploited by users without prior written consent. The information contained herein is reliable to the best of our knowledge. The indicated values are indicative or average values resulting from batches' analyses. It does not constitute any commitment nor guarantee. It is our customers' responsibility to check the local regulation before using our products.







AMINO ACIDS (g/100g protein)*					
Aspartic acid	12.3	Histidine	2.2	Proline	5.1
Glutamic acid	19.1	Isoleucine	6.0	Serine	4.8
Alanine	5.2	Leucine	13.1	Threonine	5.5
Arginine	2.8	Lysine	10.6	Tryptophane	2.5
Cystine	3.1	Methionine	2.5	Tyrosine	4.0
Glycine	2.0	Phenylalanine	4.0	Valine	5.8
*Indicative values					

STATEMENT / ADDITIONAL INFORMATION / FOOD SAFETY

GMO

Our ingredients contain no genetically modified organisms and as such are not concerned by a labelling requirement (Regulation 1829/2003/EC and 1830/2003/EC and subsequent amendments).

Various contaminants

According to our monitoring plan, Pronativ® meets the European standards or at least Codex recommendations, in terms of: heavy metals, pesticides, aflatoxins, dioxins, furans, PCBs (polychlorinated biphenyls), melamine, radioactivity...

Our powders do not contain antibiotics. They underwent neither irradiation nor ionization.

Heat treatment

According to regulation, the heat treatment of milk is at least 72 ° C/15 sec (or temperature / time equivalent).

Packaging compliance

The materials used to package our ingredients comply with European regulation EC 1935/2004 (and subsequent amendments) for materials and articles intended to be in contact with foodstuffs.

Deconditioning/Handling

We inform you that for your operations of transfer or deconditioning, and in case of accumulation of dust or powder, the Minimum Ignition Energy of our powders is between 100 and 1000 mJ and the minimum ignition or incandescence temperature (by 5 mm layer) is at least 300°C.

The explosion characteristics are given by a Kmax between 30 to 100 bar.m/s and pressure Pmax between 6.0 and 7.5 bars. Powder explosibility is between 20 and 80 g/m³. Powder resistivity is between 2.10⁹ and 3.10¹³ Ohm*m.

For particular analysis or specific regulations, please contact us.

All rights to the information contained herewith belong to Lactalis Ingredients. The information is confidential and may not be disclosed to third parties or exploited by users without prior written consent. The information contained herein is reliable to the best of our knowledge. The indicated values are indicative or average values resulting from batches' analyses. It does not constitute any commitment nor guarantee. It is our customers' responsibility to check the local regulation before using our products.